



Model GIT-P2

Model GIT-P4

Glove Integrity Tester Models

Presented by:

Dynamic Design Pharma, Inc.

Carlsbad, California USA



GIT-P2 – Dual Channel glove tester
GIT- P4 – Four channel glove tester`

GIT-P4, GIT-P2 and GIT-XA1 Glove Testers

TEST ENVIRONMENTS

SYSTEMS DESCRIPTION

OIT SCREENS, SECURITY AND DATA OUTPUT

FEATURES

PRINCIPLE OF OPERATION

INTERFACE TO ISOLATOR SYSTEM

IMPLEMENTATION NOTES

DOCUMENTATION

MODEL SELECTION

GIT-P4, GIT-P2 and GIT-XA1 Glove Testers

TEST ENVIRONMENTS

SYSTEMS DESCRIPTION

IOIT SCREENS, SECURITY AND DATA OUTPUT

FEATURES

PRINCIPLE OF OPERATION

INTERFACE TO ISOLATOR SYSTEM

IMPLEMENTATION NOTES


DOCUMENTATION

MODEL SELECTION

Negative Pressure Testing - Aseptic

Glove/Sleeve placed under vacuum from outside the isolator


In the event of a leak being present in the test item, negative pressure outside the test item causes air flow **out** of the isolator system.



Positive Pressure Testing - Containment

Glove/Sleeve is pressurized from outside the isolator

In the event of a leak being present in the test item, positive pressure inside test item causes air flow **into** the isolator system.



GIT-P4, GIT-P2 and GIT-XA1 Glove Testers

TEST ENVIRONMENTS

SYSTEMS DESCRIPTION

IOIT SCREENS, SECURITY AND DATA OUTPUT

FEATURES

PRINCIPLE OF OPERATION

INTERFACE TO ISOLATOR SYSTEM

IMPLEMENTATION NOTES


DOCUMENTATION

MODEL SELECTION



GIT-P4

- ✓ Simultaneous test of four gloves
- ✓ Stainless Steel Control Console
- ✓ Maneuvering Handle
- ✓ Polyurethane Casters
- ✓ Touchscreen Control Panel
- ✓ Printer
- ✓ Size: 500W x 630L x 1000H




GIT-P2

- ✓ Dual glove test
- ✓ Stainless Steel Control Console
- ✓ Table Top design / Optional stand
- ✓ Touchscreen Control Panel
- ✓ Printer
- ✓ Size: 400W x 400L x 250H



Control Console Rear



- ✓ Power on/off switch
- ✓ Power cord and winding bracket
- ✓ Pneumatic tubing connections
- ✓ Ethernet port
- ✓ USB connector
- ✓ ID tag



GIT-P4, GIT-P2 and GIT-XA1 Glove Testers

TEST ENVIRONMENTS

SYSTEMS DESCRIPTION

OIT SCREENS, SECURITY AND DATA OUTPUT

FEATURES

PRINCIPLE OF OPERATION

INTERFACE TO ISOLATOR SYSTEM

IMPLEMENTATION NOTES

DOCUMENTATION

MODEL SELECTION

GIT-P2 OIT Screens – Typical Navigation

GIT-P2 GLOVE INTEGRITY TESTER

CYCLE CONTROL

CH1 NOT ACTIVE

CH2 NOT ACTIVE

PARAMETERS

SYSTEM TEST

UTILITIES

DYNAMIC DESIGN PHARMA, INC.
MISSION VIEJO, CALIFORNIA - USA 05/05/14 12:43

GIT-P2 CYCLE CONTROL

	CHANNEL 1	CHANNEL 2
ENTER GLOVE ID	0	0
PRESSURE (Pa)	12345	12345
PHASE (Sec)	12345	12345
CYCLE ET (Sec)	12345	12345
STATUS	INACTIVE	INACTIVE

HOME PRINT PARAM. STOP CYCLE PRINT RESULTS

GIT-P2 GLOVE INTEGRITY TESTER

MAXIMUM PRESSURE	0 Pa
OVER PRESSURE	0 Pa
START PRESSURE	0 Pa
PASS/FAIL PRESSURE	0 Pa
INFLATE TIME	0 Sec.
STABILIZATION TIME	0 Sec.
LEAK TEST TIME	0 Sec.

SAVE ESC

GIT-P2 SYSTEM TEST





	CHANNEL 1	CHANNEL 2
PRESSURE (PV)	0 Pa	0 Pa
TIME ELAPSED	0 Sec.	0 Sec.
PRESSURE CONTROL	▲ ▼	▲ ▼

ENABLE PRINT

START TEST START TEST

HOME

GIT-P4 OIT Screens – Typical Navigation



Security Functions

- Three login levels: Operator, Supervisor and Administrator
- Individual real name assignment
- Individual Password assignment
- Automatic logout
- Login information part of the Leak Test record
- Operator has access to Leak Test functions only
- Supervisor has access to Leak test functions, Parameters Management and Utilities
- Administrator has access to all functions
- Electronic signature of each leak test record

Data Output Options

Data Output Options

- Local printing at thermal printer
- Remote printing via Ethernet
- Data storage and retrieval via laptop
- Data storage and wireless retrieval



Leak Test Data

- Active Recipe
- Time/Date Stamp
- Login user info
- Active Channel
- Parameters
- Pass/Fail Result
- Cycle Duration
- Pressure at end of test

Note:

Data output is in a .CSV file format that can be easily opened using an Excel spreadsheet

Wireless Data Transmission – All models

- DDP’s glove testers feature a wireless router that allows the password protected read-only download of the leak test data to a remote computer.
- Each day’s leak test data is automatically stored within a unique data file that is then downloaded for visualization and archiving.

19061200.CSV - Microsoft Excel

	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O	P	Q	R
1	Date	Time	OGIN leve	OGIN nam	Recipe	GLOVE ID	CYCLE #	Test Result	P (end)	ME (elapse	P (over)	P (start)	P (pass)	T (inflate)	T (stab)	(leak test	TYPE	Short Cycle
2	6/12/2019	9:04:02	ADMIN1	AD1	1	15	107	PASS	2350	334	2800	2500	2210	60	60	240	Gauntlet	OFF
3	6/12/2019	9:19:46	ADMIN1	AD1	1	15	108	FAIL	2100	344	2800	2500	2210	60	60	240	Gauntlet	OFF
4	6/12/2019	9:32:43	ADMIN1	AD1	1	15	109	PASS	2326	353	2800	2500	2210	60	60	240	Gauntlet	OFF
5	6/12/2019	9:44:48	ADMIN1	AD1	1	15	110	FAIL	2130	363	2800	2500	2210	60	60	240	Gauntlet	OFF
6	6/12/2019	9:57:45	ADMIN1	AD1	1	15	111	PASS	2338	371	2800	2500	2210	60	60	240	Gauntlet	OFF
7	6/12/2019	10:11:39	ADMIN1	AD1	1	15	112	FAIL	2130	356	2800	2500	2210	60	60	240	Gauntlet	OFF
8	6/12/2019	11:15:45	ADMIN1	AD1	1	15	113	FAIL	2161	368	2800	2500	2210	60	60	240	Gauntlet	OFF
9	6/12/2019	11:35:30	ADMIN1	AD1	1	15	114	FAIL	2196	327	2800	2500	2210	60	60	240	Gauntlet	ON
10																		

Note:

Data output is in a .CSV file format that can be easily opened using an Excel spreadsheet

Barcode Reader Functionality – All models

- DDP's glove testers feature a barcode reader that assures accuracy of the gloveport identification and eliminates the possibility of operator mistakes
- In multi-channel systems, the barcode reader also identifies the appropriate channel (A, B, C or D) by scanning the barcode located on the interface
- A “teaching” functionality allows the programming of the barcode affixed on each gloveport and each interface.



Gloveport
Scanning



Interface
Scanning

Note: Scanning of both the gloveport and the interface and automatically checking the validity of both scans assures the elimination of operator mistakes

GIT-P4 and GIT-P2 Glove Testers

TEST ENVIRONMENTS

SYSTEMS DESCRIPTION

IO SCREENS, SECURITY AND DATA OUTPUT

FEATURES

PRINCIPLE OF OPERATION

INTERFACE TO ISOLATOR SYSTEM

IMPLEMENTATION NOTES

DOCUMENTATION

MODEL SELECTION

Primary Features – All Models

- ✓ Programmable parameters
- ✓ Clear Pass/Fail test result feedback to the operator
- ✓ Numeric test result feedback
- ✓ Data Output capability
- ✓ Security login
- ✓ Positive or Negative pressure test (capability)
- ✓ Cost effective
- ✓ Barcode Reader functionality
- ✓ Simple operation, validation, training and maintenance
- ✓ Support structure for gloveport interface assemblies

Operation – All Models

STEP 1:

Operator installs the glove interface onto glove to be tested

STEP 2:

Operator starts the leak test cycle

STEP 3:

At the end of the leak test cycle, the operator reviews the results and initiates printing

STEP 4:

The system prints out the test results and the cycle concludes

NOTES:

GIT-P4 and GIT-P2 operation is identical to the GIT-XA1 operation except for the added channel selection function

GIT-P4, GIT-P2 and GIT-XA1 Glove Testers

TEST ENVIRONMENTS

SYSTEMS DESCRIPTION

IOIT SCREENS, SECURITY AND DATA OUTPUT

FEATURES

PRINCIPLE OF OPERATION

INTERFACE TO ISOLATOR SYSTEM

IMPLEMENTATION NOTES

DOCUMENTATION

MODEL SELECTION

Principle of Operation – All Models


The system detects a leak in the glove under test by comparing the internal pressure of the glove to a pass/fail pressure threshold level after having been pressurized to a set level and held in a sealed condition for a given time duration.

This methodology of leak detection is called pressure decay.


Phases of the Leak Test Cycle

- | | |
|---------------|--|
| INFLATE | Initial pressurization to set point value |
| STABILIZE | Maintain pressure at a set level for a programmed duration to allow glove material to stretch |
| DEFLATE | Allow pressure to drop to the leak test starting pressure level |
| LEAK TEST | Close off the opening to glove under test for the programmed time duration and monitor the internal pressure level. |
| LEAK TEST END | At the end of the Leak Test phase, the pressure internal to the glove under test is compared to the programmed Pass/Fail pressure level. The system then makes the leak test outcome decision and displays the result on the OIT |

Leak Test Diagram (Passing Test - Glove)



Leak Test Diagram (Failing Test / Glove)



System Capability – Glove Testing

Positive or Negative pressure, glove only

Hole size detection capability = 100 um (0.004 inches) or larger

Parameters to achieve above capability (guideline only)

- Pressure Threshold = 1500 pa
- Inflate time: less than 10 seconds
- Stabilization = 60 seconds
- Leak test = 120 seconds

Certainty of detection with > 3 sigma confidence

Certainty of no false positives with > 3 sigma confidence

System Capability – Gauntlet Testing

Positive pressure

Hole size detection capability = 150 um (0.006 inches) or larger


Cycle time parameters to achieve above capability (guideline only)

- Inflate time: 60 seconds (approx.)
- Stabilization = 120 seconds
- Leak test = 180 seconds

Certainty of detection with > 3 sigma confidence

Certainty of no false positives with > 3 sigma confidence

Process window determination



GIT-P4, GIT-P2 and GIT-XA1 Glove Testers

TEST ENVIRONMENTS

SYSTEMS DESCRIPTION

IO SCREENS, SECURITY AND DATA OUTPUT

FEATURES

PRINCIPLE OF OPERATION

INTERFACE TO ISOLATOR SYSTEM

IMPLEMENTATION NOTES

DOCUMENTATION

MODEL SELECTION

Interface to the Isolator - Options



Containment Applications (Positive Pressure Testing)

1. Glove Interface
2. Single Piece Gauntlet Interface


Aseptic Applications (Negative Pressure Testing)

1. Glove Vacuum Chamber
2. Single Piece Gauntlet Vacuum Chamber


Glove Interface - Positive Pressure Testing





Gauntlet Interface - Positive Pressure Testing



Glove Vacuum Chamber - Negative Pressure Testing



Gauntlet Vacuum Chamber – Negative Pressure Testing



GIT-P4, GIT-P2 and GIT-XA1 Glove Testers

TEST ENVIRONMENTS

SYSTEMS DESCRIPTION

IOIT SCREENS, SECURITY AND DATA OUTPUT

FEATURES

PRINCIPLE OF OPERATION

INTERFACE TO ISOLATOR SYSTEM

IMPLEMENTATION NOTES

DOCUMENTATION

MODEL SELECTION

Validation Principle

- Obtain a glove known to be leak free (multiple gloves can be used as well)
- Run a minimum of 10 separate tests with and without the DDP supplied test orifice connected
- Determine whether an acceptable process window is present

GIT-P4, GIT-P2 and GIT-XA1 Glove Testers

TEST ENVIRONMENTS

SYSTEMS DESCRIPTION

IO SCREENS, SECURITY AND DATA OUTPUT

FEATURES

PRINCIPLE OF OPERATION

INTERFACE TO ISOLATOR SYSTEM

IMPLEMENTATION NOTES

DOCUMENTATION

MODEL SELECTION

Documentation

The following documentation package is supplied with the system

- Instruction manual (operator and maintenance)
- Assembly and control drawing package
- System Specification
- Functional and Design Specification document
- Functional Test document (executed)
- Factory Acceptance Test (executed)
- IOQ Protocol (ready for execution)

GIT-P4, GIT-P2 and GIT-XA1 Glove Testers

TEST ENVIRONMENTS

SYSTEMS DESCRIPTION

IOIT SCREENS, SECURITY AND DATA OUTPUT

FEATURES

PRINCIPLE OF OPERATION

INTERFACE TO ISOLATOR SYSTEM

IMPLEMENTATION NOTES

DOCUMENTATION

MODEL SELECTION

GIT-P4, GIT-P2 and GIT-S1 Glove Testers

	GIT-P4	GIT-P2	GIT-S1
Number of gloves simultaneously tested	4	2	1
Test methodology	Press. Decay	Press. Decay	Press. Decay
Number of gloves tested/hour (including set up)	80 (3)(5)	40 (3)(5)	15 (1)
Number of gauntlets tested/hour (including set up)	40 (4)(5)	20 (4)(5)	8 (1)
Recommended isolator size (# of gloves)	>12	6-12	2-6

- (1) This calculation assumes a gloveport interface installation time of 1 minute.
- (2) All models require a single electrical power connection
- (3) The above assumes a glove testing cycle time of 3 minutes
- (4) The above assumes a gauntlet testing cycle time of 6 minutes
- (5) This calculation assumes no gloveport interface installation time (continuous testing)
- (6) The GIT-S1 tester is not described in this presentation

Summary

- Clean room friendly, easy to operate and maintain
- Capable of meeting hole size detection requirements
- Can be validated using non-microbiological methods
- Highly cost effective

Thank You

Dynamic Design Pharma